


ISOLATE YOUR EYES  
**COCOONS®**

**CLIP-ONS SIZING INSTRUCTIONS:**

Use the shape illustration to determine if this clip size most accurately fits your prescription glasses. To use the sizing guide, print out this form at 100% actual size and then place your prescription glasses face down on the shape to determine if it fits your eyewear.

This clip fits a maximum of 120 mm (4 3/4 inches).

42 mm  
(1 5/8 inches)



● RND2-48

This clip fits a maximum of 126 mm (4 15/16 inches).

44 mm  
(1 3/4 inches)



● RND2-50

This clip fits a maximum of 134 mm (5 1/4 inches).

46 mm  
(1 3/4 inches)



● RND2-52

This clip fits a maximum of 139 mm (5 1/2 inches).

48 mm  
(1 7/8 inches)



● RND2-54



INCLUDED CASE



The image of the quarter should be the same size as an actual quarter when printed at full scale. The image should be 24 mm in diameter at full size.